

Diversey™

# Endbac Sanitising wipes

## Sanitising wipes

- Broad spectrum bactericide
- Maximised bacteria kill
- Passes EN1276
- Efficient and convenient

### Maximised bacteria kill

Endbac probe wipes are effective across a wide range of common food poisoning bacteria including E.Coli and S. Aureus. The product formula incorporates a preservative system that ensures effective germ kill is maintained for up to two years if the product is stored correctly.

### Efficient and convenient

Endbac probe wipes are made from a high quality, non-woven material that offers excellent strength and absorbency. Because they are designed for single use and are disposable, they are ideal for “clean as you go” tasks.

Endbac probe wipes are suitable for sanitising probes and utensils between applications.





# Endbac Sanitising wipes

The professional choice for kitchen sanitising

## Usage Instructions

- Pull each wipe sharply at an angle to dispense.
- Disinfect the probe with a careful, firm wiping action.
- This product is suitable for use on keyboards, telephones, desks and other hard surfaces
- Also suitable for disinfecting utensils.
- Wipe surface with wet wipe and leave it to dry, rinse the (cleansed) surface / utensil before use.
- After use, dispose of each wipe - do not re-use.

## Technical data

Appearance	Blue woven wipes	pH	7.0
------------	------------------	----	-----

*The above data is typical of normal production and should not be taken as a specification.*

## Packaging Information

Product Code	Pack Size	Packs per case	Cases per layer	Cases per pallet	Outer Bar Code	Inner Bar Code
101100971	200pc	6	18	108	7615400778492	7615400778508
101100972	1000pc		20	80	7615400778515	